

medax

MEDWIRE & SINUREP

WIRE LOCALIZATION NEEDLES FOR THE BREAST



medax
medical devices

Your partner in Biopsy and Special Needles

MEDAX


MEDWIRE

REPOSITIONABLE WIRE LOCALIZATION NEEDLE FOR NON-PALPABLE BREAST LESIONS

- Easy “button” wire release
- Extremely sharp tip for atraumatic penetration
- Echogenic treatment for ultrasound guided procedures
- Centimeter markings and sliding stopper
- Depth indicator on the wire
- 5 year shelf life



O-TYPE HOOK WIRE



J-TYPE HOOK WIRE

	GAUGE Ø	WIRE (MM)	LENGTH (MM)		
O-Type Hook Wire	20	0,40	20 060	20 090	20 120
J-Type Hook Wire	20	0,40	20 030	20 050	20 070 20 100

Other lengths available as special orders; minimum order quantities may apply

Ordering information:

MEDWIRE O-TYPE HOOK WIRE

M	W	GAUGE	LENGTH	2	0
---	---	-------	--------	---	---

MEDWIRE J-TYPE HOOK WIRE

M	W	GAUGE	LENGTH	1	1
---	---	-------	--------	---	---

Single box: 10 PCS.



K.M. Medical, Inc.

65 Parker St. Unit 2, Newburyport, MA 01950

P (978) 465-4689 • F (978) 914-6106

info@kmmedical.com


Excellence through Customer Satisfaction

medax

SINUREP

WIRE LOCALIZATION NEEDLE FOR NON-PALPABLE BREAST LESIONS

- Extremely sharp tip for atraumatic penetration
- Echogenic treatment for ultrasound guided procedures
- Centimeter markings and sliding stopper
- 5 year shelf life



VERSION HOOK	GAUGE Ø	WIRE (MM)	LENGTH (MM)		
L-Single	20	0,25-0,30	20 080	20 100	20 120

Other lengths and wire configuration available as special orders; minimum quantities may apply

Ordering information:

SINUREP SINGLE HOOK WIRE

S	R	GAUGE	LENGTH	0	0
---	---	-------	--------	---	---

SINUREP SINGLE HOOK WIRE MARKED

S	R	GAUGE	LENGTH	0	2
---	---	-------	--------	---	---

SINUREP SINGLE HOOK WIRE STIFFENED VERSION

S	R	GAUGE	LENGTH	0	3
---	---	-------	--------	---	---



Excellence through Customer Satisfaction

K.M. Medical, Inc.

65 Parker St. Unit 2, Newburyport, MA 01950

P (978) 465-4689 ▪ F (978) 914-6106

info@kmmedical.com



Your partner in Biopsy and Special Needles

For ordering information, contact KM Medical Sales specialist or Customer Service.



Excellence through Customer Satisfaction

K.M. Medical, Inc.
65 Parker St. Unit 2,
Newburyport, MA 01950
P (978)465-4689 ▪ F (978)914-6106
info@kmmedical.com



Medax's quality management system is certified to ISO 13485:2012 standards, in accordance with the requirements of the medical device 93/42/CEE directive and its relevant updates. All products undergo intensive clinical testing and are fully CE and FDA approved.