

Tumark[®]



Marker System




Excellence through Customer Satisfaction

K.M. Medical, Inc.
65 Parker St. Unit 2
Newburyport, MA 01950
P: 978-465-4689 F: 978-914-6106
Email: info@kmmmedical.com

Tumark® Q

The Tumark® Q serves for percutaneous marking in the pectoral tissue. Among other things, the marker is implanted in lesions prior to or during chemotherapy and the location of a biopsy sampling point. The marker is visible under mammography, ultrasound and MRI *

* Only the marker of Tumark® Q is MRI compatible.



Features:

- Q-merged form enables firm anchorage in the breast tissue and increases visibility in all positions
- Visibility of the marker in ultrasound, mammography and MRI
- Marker is MRI compatible up to 3 Tesla
- Marker is approved as implantable material (Nitinol)
- Marker is supplied preloaded
- Ergonomic handle for one single-hand operation
- Simultaneous use with Tumark Professional (U-Marker) enables better differentiation during mammography

Tumark® Q

Ref. no.	Gauge	Marker Size	Diameter	Length
271 500	18G	4,2 x 4 mm	1,2 mm	120 mm

Fully MRI-compatible

DS-5-11 Version 4



Excellence through Customer Satisfaction

K.M. Medical, Inc.
 65 Parker St. Unit 2
 Newburyport, MA 01950
P: 978-465-4689 **F:** 978-914-6106
Email: info@knmedical.com