



## FEATURES

- 11 different sizes ranging from 10x5mm to 25x10mm
- Loop, ball, and square electrodes

## TECHNICAL DATA

Shaft Length . . . . . 13cm (5")  
Sterilization Method . . . . . ETO  
Box Qty. . . . . 20 ea.

## ORDERING INFORMATION

### SQUARE ELECTRODES

LE-06-108 - 10 x 8mm, Square Electrode  
LE-06-101 - 10 x 10mm, Square Electrode

### BALL ELECTRODES

LE-13-003 - 3mm Diameter, Ball Electrode  
LE-14-005 - 5mm Diameter, Ball Electrode

### LOOP ELECTRODES

LE-04-105 - 10 x 5mm, Loop Electrode  
LE-05-101 - 10 x 10mm, Loop Electrode  
LE-07-158 - 15 x 8mm, Loop Electrode  
LE-26-150 - 15 x 10mm, Loop Electrode  
LE-09-208 - 20 x 8mm, Loop Electrode  
LE-24-201 - 20 x 10mm, Loop Electrode  
LE-11-205 - 20 x 15mm, Loop Electrode



Download IFU



Made in the USA  
ISO 13485 and FDA Registered



65 Parker St., Unit #2,  
Newburyport, MA 01950  
Phone: (800) 531-2545

Fax: (978) 914-6106  
Email: [info@kmmedical.com](mailto:info@kmmedical.com)