

**medax**

# MEDEASY

SEMI-AUTOMATIC SPRING LOADED BIOPSY SYSTEM




**K.M. Medical, Inc.**  
65 Parker St. Unit 2, Newburyport, MA 01950  
P (978) 465-4689 ■ F (978) 914-6106  
info@kmmedical.com

*Excellence through Customer Satisfaction*

**MEDAX**

**MEDEASY**

SEMI-AUTOMATIC SPRING LOADED BIOPSY SYSTEM



Semi-automatic spring loaded biopsy system with selectable penetration depth.

- Adjustable penetration depth (10/20 mm).
- Echogenic treatment for ultrasound guided procedures.
- Centimeter markings and optional sliding stopper.
- Also available with Sword tip.
- Can be used with CE introducer coaxial needle.
- Sterilized by ETO, shelf life 5 years.



Medax's quality management system is certified to ISO 13485: 2016 standards, in accordance with the requirements of the medical device 93/42/CEE directive and its relevant updates. All products undergo intensive clinical testing and are fully CE and FDA approved.



**K.M. Medical, Inc.**  
65 Parker St. Unit 2, Newburyport, MA 01950  
P (978) 465-4689 ■ F (978) 914-6106  
info@kmmmedical.com

*Excellence through Customer Satisfaction*

**MEDAX**

# MEDEASY

## SEMI-AUTOMATIC SPRING LOADED BIOPSY SYSTEM

GAUGE	COLOUR	LENGTH (MM)							
		12 080	12 100	12 120	12 160	12 180	12 200	12 250	12 300
12	● Light blue	12 080	12 100	12 120	12 160	12 180	12 200	12 250	12 300
14	● Green	14 080	14 100	14 120	14 160	14 180	14 200	14 250	14 300
16	○ White	16 080	16 100	16 120	16 160	16 180	16 200	16 250	16 300
18	● Pink	18 080	18 100	18 120	18 160	18 180	18 200	18 250	18 300
20	● Yellow	20 080	20 100	20 120	20 160	20 180	20 200	20 250	20 300

Different sizes on demand

### Ordering information:

MEDEASY BEVEL TIP

M	S	GAUGE		LENGTH		0	0

MEDEASY BEVEL TIP "LIGHT VERSION"

M	L	GAUGE		LENGTH		0	0

MEDEASY SWORD TIP

M	S	GAUGE		LENGTH		0	2

MEDEASY SWORD TIP "LIGHT VERSION"

M	L	GAUGE		LENGTH		0	2

INTRODUCER COAXIAL NEEDLE FOR MEDEASY

C	E	GAUGE		LENGTH		0	0

MEDEASY BEVEL TIP WITH COAXIAL NEEDLE

M	S	GAUGE		LENGTH		C	0

MEDEASY BEVEL TIP "LIGHT VERSION" WITH COAXIAL NEEDLE

M	L	GAUGE		LENGTH		C	0

MEDEASY BEVEL TIP WITH COAXIAL NEEDLE INCLUDING BLUNT TIP STYLET

M	S	GAUGE		LENGTH		K	0

MEDEASY BEVEL TIP "LIGHT VERSION" WITH COAXIAL NEEDLE INCLUDING BLUNT TIP STYLET

M	L	GAUGE		LENGTH		K	0

Single box: 10 PCS.



**K.M. Medical, Inc.**  
 65 Parker St. Unit 2, Newburyport, MA 01950  
 P (978) 465-4689 ■ F (978) 914-6106  
 info@kmmmedical.com

*Excellence through Customer Satisfaction*



Your partner in Biopsy and Special Needles

For ordering information, contact KM Medical Sales specialist or Customer Service.



*Excellence through Customer Satisfaction*

K.M. Medical, Inc.  
65 Parker St. Unit 2,  
Newburyport, MA 01950  
P (978)465-4689 ▪ F (978)914-6106  
info@kmmedical.com



Medax's quality management system is certified to ISO 13485:2012 standards, in accordance with the requirements of the medical device 93/42/CEE directive and its relevant updates. All products undergo intensive clinical testing and are fully CE and FDA approved.