

MEDaX

MEDONE

AUTOMATIC DISPOSABLE BIOPSY SYSTEM



- One hand usage
- 22mm penetration depth, 19mm sample notch
- Dual firing triggers in proximal and distal positions
- Echogenic tip for ultrasound guided procedures
- Centimeter markings
- Can be used with coaxial introducer needle
- 5 year shelf life

| GAUGE | COLOR CODE | LENGTH (MM) | | | |
|-------|--------------|-------------|--------|--------|---------------|
| 12 | ● Light blue | 12 100 | 12 160 | | |
| 14 | ● Green | 14 100 | 14 160 | 14 200 | |
| 16 | ○ White | 16 100 | 16 160 | 16 200 | |
| 18 | ● Pink | 18 100 | 18 160 | 18 200 | 18 250 18 300 |
| 20 | ● Yellow | 20 100 | 20 160 | 20 200 | |

Other lengths and tips available as special orders; minimum quantities may apply

Ordering information:

MEDONE BEVEL TIP

| | | | | | |
|---|---|-------|--------|--|--|
| M | O | GAUGE | LENGTH | | |
| | | | | | |

COAXIAL INTRODUCER NEEDLE FOR MEDONE

| | | | | | |
|---|---|-------|--------|--|--|
| C | O | GAUGE | LENGTH | | |
| | | | | | |

MEDONE WITH COAXIAL NEEDLE IN THE SAME PACKAGE

| | | | | | |
|---|---|-------|--------|---|---|
| M | O | GAUGE | LENGTH | C | O |
| | | | | | |

Single Box: 10 PCS



Excellence through Customer Satisfaction

K.M. Medical, Inc.
65 Parker St. Unit 2, Newburyport, MA 01950
P (978) 465-4689 • F (978) 914-6106
info@kmmedical.com

medaX

MEDONE

AUTOMATIC DISPOSABLE BIOPSY SYSTEM

| GAUGE | COLOUR | LENGTH (MM) | | | | | | | |
|-------|--------------|-------------|--------|--------|--------|--------|--------|--------|--------|
| | | 080 | 100 | 120 | 160 | 180 | 200 | 250 | 300 |
| 12 | ● Light blue | 12 080 | 12 100 | 12 120 | 12 160 | 12 180 | 12 200 | 12 250 | 12 300 |
| 14 | ● Green | 14 080 | 14 100 | 14 120 | 14 160 | 14 180 | 14 200 | 14 250 | 14 300 |
| 16 | ○ White | 16 080 | 16 100 | 16 120 | 16 160 | 16 180 | 16 200 | 16 250 | 16 300 |
| 18 | ● Pink | 18 080 | 18 100 | 18 120 | 18 160 | 18 180 | 18 200 | 18 250 | 18 300 |
| 20 | ● Yellow | 20 080 | 20 100 | 20 120 | 20 160 | 20 180 | 20 200 | 20 250 | 20 300 |

Different sizes on demand

Ordering information:

MEDONE NORMAL TIP

| | | | | | |
|---|---|-------|--------|---|---|
| M | O | GAUGE | LENGTH | 0 | 0 |
| | | | | | |

MEDONE TROCAR TIP

| | | | | | |
|---|---|-------|--------|---|---|
| M | O | GAUGE | LENGTH | 0 | 1 |
| | | | | | |

INTRODUCER COAXIAL NEEDLE FOR MEDONE

| | | | | | |
|---|---|-------|--------|---|---|
| C | O | GAUGE | LENGTH | 0 | 0 |
| | | | | | |

MEDONE WITH COAXIAL NEEDLE IN THE SAME PACKAGE

| | | | | | |
|---|---|-------|--------|---|---|
| M | O | GAUGE | LENGTH | C | 0 |
| | | | | | |

MEDONE WITH COAXIAL NEEDLE INCLUDING BLUNT TIP STYLET
IN THE SAME PACKAGE

| | | | | | |
|---|---|-------|--------|---|---|
| M | O | GAUGE | LENGTH | K | 0 |
| | | | | | |

Single box: 10 PCS.

*Excellence through Customer Satisfaction***K.M. Medical, Inc.**

65 Parker St. Unit 2, Newburyport, MA 01950

P (978) 465-4689 ▪ F (978) 914-6106

info@kmmmedical.com



Your partner in Biopsy and Special Needles

For ordering information, contact KM Medical Sales specialist or Customer Service.



Excellence through Customer Satisfaction

K.M. Medical, Inc.
65 Parker St. Unit 2,
Newburyport, MA 01950
P (978)465-4689 ▪ F (978)914-6106
info@kmmedical.com



Medax's quality management system is certified to ISO 13485:2012 standards, in accordance with the requirements of the medical device 93/42/CEE directive and its relevant updates. All products undergo intensive clinical testing and are fully CE and FDA approved.