

medax

MEDCUT

AUTOMATIC DISPOSABLE BIOPSY SYSTEM



Automatic disposable spring loaded biopsy system.

- 11/22 mm penetration depth.
- Firing trigger in proximal position.
- Sharpened cannula for a safe and painless percutaneous insertion and penetration.
- Echogenic treatment for ultrasound guided procedures.
- Centimeter markings and optional sliding stopper.
- Also available with the Trocar tip.
- Can be used with CM introducer coaxial needle.
- Sterilized by ETO, shelf life 5 years.



Excellence through Customer Satisfaction

K.M. Medical, Inc.
65 Parker St. Unit 2, Newburyport, MA 01950
P (978) 465-4689 • F (978) 914-6106
info@kmmedical.com

MEDCUT

AUTOMATIC DISPOSABLE BIOPSY SYSTEM

GAUGE	COLOUR	LENGTH (MM)							
12	● Light blue	12 080	12 100	12 120	12 160	12 180	12 200	12 250	12 300
14	● Green	14 080	14 100	14 120	14 160	14 180	14 200	14 250	14 300
16	○ White	16 080	16 100	16 120	16 160	16 180	16 200	16 250	16 300
18	● Pink	18 080	18 100	18 120	18 160	18 180	18 200	18 250	18 300
20	● Yellow	20 080	20 100	20 120	20 160	20 180	20 200	20 250	20 300

Different sizes on demand

Ordering information:

MEDCUT NORMAL TIP (penetration depth 22 mm) M C GAUGE LENGTH 0 0

MEDCUT TROCAR TIP (penetration depth 22 mm) M C GAUGE LENGTH 0 1

MEDCUT SHORT NOTCH NORMAL TIP (penetration depth 11 mm) M C GAUGE LENGTH 0 S

INTRODUCER COAXIAL NEEDLE FOR MEDCUT (penetration depth 22 mm) C M GAUGE LENGTH 0 0

INTRODUCER COAXIAL NEEDLE FOR MEDCUT (penetration depth 11 mm) C M GAUGE LENGTH 0 S

MEDCUT NORMAL TIP (penetration depth 22 mm)
WITH COAXIAL NEEDLE IN THE SAME PACKAGE M C GAUGE LENGTH C 0

MEDCUT SHORT NOCH NORMAL TIP (penetration depth 11 mm)
WITH COAXIAL NEEDLE IN THE SAME PACKAGE M C GAUGE LENGTH C S

Single box: 10 PCS.





Your partner in Biopsy and Special Needles

For ordering information, contact KM Medical Sales specialist or Customer Service.



Excellence through Customer Satisfaction

K.M. Medical, Inc.
65 Parker St. Unit 2,
Newburyport, MA 01950
P (978)465-4689 ▪ F (978)914-6106
info@kmmedical.com



Medax's quality management system is certified to ISO 13485:2012 standards, in accordance with the requirements of the medical device 93/42/CEE directive and its relevant updates. All products undergo intensive clinical testing and are fully CE and FDA approved.